

PROTOCOL EXCHANGE | COMMUNITY CONTRIBUTED Immunofluorescence of Cultured Cells

Paige Taylor

Krakow Lab

Abstract

This protocol describes the PFA-fixation and immunofluorescent-labeling of cultured cells.

Subject terms: Cell biology Cell culture Imaging

Keywords: immunofluorescence

Reagents

4% Paraformaldehyde

PBS

PBS + 0.1% TritonX-100

10% serum (from species of secondary)

Primary antibody

Alexa Fluor 488 secondary

Alexa Fluor 568 secondary

DAPI

Mounting media

Equipment

Humidity Chamber

Procedure

1. Wash cells with PBS.
2. Fix in 4% PFA for 15 min at room temperature (RT).
3. Wash cells with PBS (2×5 min). Cells can be stored in PBS at 4C for up to a week.
4. Permeabilize with PBS + 0.1% TritonX-100 for 5 min at RT.
5. Wash with PBS (2×5 min) at RT.
6. Block with 10% serum (from species of secondary) for 1 hr at RT.
7. Drain and incubate with primary antibody diluted in 5% serum overnight at 4C.
8. Wash in PBS (3×5 min) at RT.

9. Detect primary antibody with Alexa Fluor secondary antibody diluted 1:400 in 5% serum at RT for 1 hr in the dark.
10. Wash in PBS (3×5 min) at RT.
11. Add mounting media with DAPI and apply coverslip.
12. Store slides in the dark at 4C.

Author information

Affiliations

Krakow Lab

Paige Taylor

Competing financial interests

None.

Corresponding author

Correspondence to: Paige Taylor (sptaylor@g.ucla.edu)

Readers' Comments

Comments on this thread are vetted after posting.

Protocol Exchange ISSN 2043-0116

© 2015 Nature Publishing Group, a division of Macmillan Publishers Limited. All Rights Reserved.
partner of AGORA, HINARI, OARE, INASP, CrossRef and COUNTER