

PROTOCOL EXCHANGE | COMMUNITY CONTRIBUTED ^{ED} Reprogramming mouse embryonic fibroblasts using different reprogramming factors

Jing Liu , Jiekai chen & Duanqing Pei

Pei's Lab,GIBH

Abstract

MEFs can be reprogramed to pluripotent stem cells by defined reprogramming factors¹. However, the reprograming efficiency is influenced by kinds of conditions, such as the combination of reprogramming factors, culture medium, types of the virus and the cell density etc². This protocol described the optimized condition for reprogramming MEFs by different reprogramming system, especially the combination of different reprogramming factors.

Subject terms: Cell biology Cell culture

Keywords: Somatic reprogramming mouse embryonic fibroblast
Reprogramming factors

Reagents

- ☐ DMEM high glucose (HyClone; #SH30022.01B)
- ☐ FBS(MEF culture) (PAA; #A15-101)
- ☐ FBS (ES qualify) (GIBCO; #10099141)
- ☐ Glutamax (GIBCO ;#35050079)
- ☐ non-essential amino acids (GIBCO; #11140076)
- ☐ sodium pyruvate (GIBCO ;#11360070)
- ☐ β -mercaptoethanol(GIBCO;#21985023)
- ☐ LIF (Enzo Life Science ;#ALX-201-242)
- ☐ Ascorbic acid(sigma ;#A4403)
- ☐ Gelatin(sigma;#119k0062)
- ☐ Trypsin-EDTA (Life Technologies; 25300054)
- ☐ Recombinant human fibroblast growth factor (FGF2) (Peprotech; 100-18B)
- ☐ Recombinant human epidermal growth factor (EGF) (Peprotech; AF-100-15)
- ☐ Recombinant human BMP4(R&D;#314-Bp)
- ☐ DMSO(Sigma;SHBC4445V)

- ☐ Y27632 (Calbiochem; 688001)
- ☐ TTNPB(Selleck; S4627)
- ☐ DZNep(Selleck; S7120)
- ☐ RepSox(Sigma; R0158)
- ☐ CHIR99021(MCE; HY-10182)
- ☐ PD0325901(MCE;HY-10254)
- ☐ VPA – Valproic acid sodium salt (Sigma; P4543)

Procedure

pMXs-Based Retrovirus induced reprogramming

Preparation of the Virus

- 1) Plate 8×10^6 Plat-E cells per 10-cm dish uniformly. Cells were cultured in 10%FBS medium for 16 hours to reach a 70-80% confluent.
- 2) Transfection was carried out with the modified calcium phosphate transfection method³ as follows: Replace the medium of the Plan-E cells in 10-cm dish with 7.5 ml fresh 10% FBS medium. For each factor, 1068 μ l ddH₂O, 25 μ g plasmid, 156.25 μ l 2M CaCl₂, 1.25ml 2XHBS(total 2.5ml) were added to a 15ml tube successively. Mix vigorously after adding the 2XHBS, and then incubate for 2 min at Room temperature. Transfer the 2.5 ml mixture into the dish that had replaced with 7.5ml fresh medium above.
- 3) 11-14 hours after transfection, replace the medium of the transfected 10 cm plat-E cell dishes with 10ml fresh medium. Continue to incubate the cells.
- 4) 48 hours after transfection, the Supernatant contained the virus was collected by a syringe and filter though a 0.45 μ m filter as the first virus stock. A 10ml fresh Plan-E medium was added to the transfected Plat-E cells and harvest 24 hours later as the second virus stock.

Preparation of the OG2 MEFs

Thawing the frozen Passage 1 OG2 MEF into a 6 cm dish with MEF Medium, and cultured in the CO₂ incubator to reach a 100% confluence. Then split the MEFs to P12 or P24 plate at different cell density.

NOTE: As the cell plated density have a significant influence to the reprogramming efficiency². We adjusted the plated density according to the combination of reprogramming factors. In summary, a higher density was used for a lower reprogramming efficiency system and a higher density was used for a lower reprogramming system. For detail:

- ☐ O/K/S/M: 1.2-1.5X10⁴ cells per well of 12-well plate
- ☐ O/K/S:1.5-2 X10⁴ cells per well of 12-well plate
- ☐ KSM/OSM/OKM: 1.5-2 X10⁴ cells per well of 12-well plate
- ☐ c-JunDN/K/S: 2 X10⁴ cells per well of 12-well plate
- ☐ Jdp2/K/S: 2 X10⁴ cells per well of 12-well plate

- ☐ O/SI/d1: 3 X104 cells per well of 12-well plate
- ☐ Oct4/Id1/Jhdm1b/Lrh1/Glis1/Sall4: 3 X104 cells per well of 12-well plate

Infect the MEFs with the virus stock

1) Mix the virus stock with equal volume and add one volume fresh MEF medium.

For example:

OKSM induced reprogramming: 0.5ml Oct4, 0.5ml Klf4, 0.5ml Sox2, 0.5ml Myc, and 0.5ml fresh MEF medium (Total 2.5ml) for one well MEFs of 12-well plate.

OKS induced reprogramming: 0.5ml Oct4, 0.5ml Klf4, 0.5ml Sox2, and 0.5ml fresh MEF medium (Total 2ml) for one well MEFs of 12-well plate.

OKS induced reprogramming: 0.25ml Oct4, 0.25ml Klf4, 0.25ml Sox2, and 0.25ml fresh MEF medium (Total 2ml) for one well MEFs of 24-well plate.

Oct4/Jhdm1b/Id1 induced reprogramming: 0.5ml Oct4, 0.5ml Jhdm1b, 0.5ml Id1, and 0.5ml fresh MEF medium (Total 2ml) for one well MEFs of 12-well plate.

Jdp2/Id1/Jhdm1b/Lrh1/Glis1/Sall4: 0.5ml JDP2, 0.5ml Jhdm1b, 0.5ml Id1, 0.5ml Lrh1, 0.5ml Id1, 0.5ml Sall4 and 0.5ml fresh MEF medium (Total 3.5ml) for one well MEFs of 12-well plate.

2) Add polybrene to a final concentration of 4ug/ml.

3) Repeat the infection with the second virus stock 24 hours after the first infection.

The Day after the second infection noted as post-infection Day0.

Generation of the iPSCs

1) At post-infection Day0, replace the virus contained medium with fresh reprogramming medium such as mES , mES+Vc4 or iCD15 medium., according to the design of the experiment. (1ml per 12-plate well)

2) Change the medium everyday and observe the morphology change.

3) Count the GFP+ colonies at proper days.

4) Pick the GFP+ colonies and passage.

pSuper bashed shRNA induced reprogramming

All of the protocols are same with the pMXs-Based Retrovirus except that, the puromycin (2µg/ml) was added to the reprogramming medium at post-infection Day2 to root out the uninfected cells.

pW-TRE based Doxycycline(Dox) induced reprogramming.

Preparation of the Virus

1) Plate 6×10⁶ HEK93T cells per 10-cm dish uniformly. Cells were cultured in 10%FBS medium for 16 hours to reach a 70-80% confluent.

2) Transfection was carried out with the modified calcium phosphate transfection method as follows: Replace the medium of the HEK93T cells in 10-cm dish with 7.5 ml fresh 10% FBS medium. For each factor, 1070µl ddH₂O, 12.5µg pW-TRE-c-Jun/ pW-rtTA, 7.5µg psPAX2, 5µg

psMD2.G, 156.25µl 2M CaCl₂, and 1.25ml 2XHBS(total 2.5ml), were added to a 15ml tube successively. Mix vigorously after adding the 2XHBS, and then incubate for 2 min at Room temperature. Transfer the 2.5 ml mixture into the dish that had replaced with 7.5ml fresh medium above.

3) 11-14 hours after transfection, replace the medium of the transfected 10 cm HEK93T cell dishes with 10ml fresh medium. Continue to incubate the cells.

4) 48 hours after transfection, the Supernatant contained the virus was collected by a syringe and filter through a 0.45 µm filter as the first virus stock. A 10ml fresh HEK93T medium was added to the transfected HEK93T cells and harvest 24 hours later as the second virus stock.

Preparation of the MEFs

Same with relative section of “pMXs-Based Retrovirus induced reprogramming”

Infect the MEFs with virus stock

Infect the MEFs with virus carrying the Oct4, Klf4, Sox2 and c-Jun coding sequences.

To reprogram cells in a P12 well:

1) mix the virus stock as following:

0.5ml Oct4 (pMXs-Based Retrovirus), 0.5ml Sox2 (pMXs-Based Retrovirus), 0.5ml Klf4 (pMXs-Based Retrovirus), 0.125ml c-Jun (pW-TRE-Based Lentivirus), 0.125ml rtTA (pW-Based Lentivirus), and 0.5ml fresh MEF medium.

2) Add polybrene to a final concentration of 4ug/ml.

3) Repeat the infection with the second virus stock 24 hours after the first infection.

The Day after the second infection noted as post-infection Day0.

Generation of the iPSCs

1) At post-infection Day0, replace the virus contained medium with fresh reprogramming medium such as mES, mES+Vc or iCD1 medium. DOX was added in the indicated days as design of the experiment.

2) Change the medium everyday and observe the morphology change.

3) Count the GFP+ colonies at proper days.

Timing

Virus stock preparation: 3 Days

MEFs preparation: 3 Days

OKS induced reprogramming in iCD1 medium: 8 Days

OKS/c-JunDN induced reprogramming in iCD1 medium: 8 Days

OKS/c-Jun shRNA induced reprogramming in iCD1 medium: 10 Days

OKS induced reprogramming in mES medium: 21days

OKS induced reprogramming in mES+Vc medium: 16 days

OKSM induced reprogramming in iCD1 medium: 6 Days

OKSM induced reprogramming in mES+Vc medium: 12 Days

c-JunDN/KSM, JDP2/KSM induced reprogramming in iCD1 medium medium: 8 Days

c-Jun shRNA/KSM induced reprogramming in iCD1 medium medium: 8 Days

c-JunDN/KS, JDP2/KS induced reprogramming in iCD1 medium medium: 16 Days

Oct4/Sox induced reprogramming in iCD1+BMP medium: 16 Days.

Oct4/Jhdm1b/Id1 induced reprogramming in iCD1+BMP medium: 14 Days.

Jdp2/Id1/Jhdm1b/Lrh1/Glis1/Sall4 induced reprogramming in iCD1 medium: 16 Days.

Troubleshooting

Problem 1 Cells grow slowly.

Reason: 1)MEFs are over grow before split; 2)The FBS is not qualified.

Solution: 1)Use MEFs in passage number 2 and split the MEFs once the confluence reach to 100% after thawing; 2)Use the qualified FBS.

Problem 2 The morphology change of the reprogramming cell is slowly

Reason: 1) The titer of the virus is low; 2)The medium is out of date.

Solution: 1)repare the virus following the protocol strictly; 2) Use fresh virus supernatant; 3) Use qualified reprogramming medium.

Anticipated Results

For OKS induced reprogramming in iCD1 medium, a significant mesenchyme to epithelial transition can be observed at post-infection Day2, and GFP+ cells will appear at post-infection Day4. More than 800 GFP+ colonies can be found at post-infection Day8 in a P12 well. c-Jun over expression can significantly inhibit the mesenchyme to epithelial transition induced by OKS, and few colonies can be found at post-infection Day8. For Jdp2/Id1/Jhdm1b/Lrh1/Glis1/Sall4 induced reprogramming in iCD1 medium, cells grow at a relative low rate, and mesenchyme to epithelial transition can be found at post-infection Day6, GFP+ colonies can be found at post-infection Day16.

References

1. Takahashi, K. & Yamanaka, S. Induction of pluripotent stem cells from mouse embryonic and adult fibroblast cultures by defined factors. *Cell* 126, 663-676 (2006).
2. Chen, J. et al. Towards an optimized culture medium for the generation of mouse induced pluripotent stem cells. *The Journal of biological chemistry* 285, 31066-31072 (2010).
3. Urabe, M., Kume, A., Tobita, K. & Ozawa, K. DNA/Calcium phosphate precipitates mixed with medium are stable and maintain high transfection efficiency. *Analytical biochemistry* 278, 91-92 (2000).

4. Esteban, M.A. et al. Vitamin C enhances the generation of mouse and human induced pluripotent stem cells. *Cell stem cell* 6, 71-79 (2010).
5. Chen, J. et al. Rational optimization of reprogramming culture conditions for the generation of induced pluripotent stem cells with ultra-high efficiency and fast kinetics. *Cell research* 21, 884-894 (2011).

Acknowledgements

We thanks ail members of Pei labs for their support.

Associated Publications

This protocol is related to the following articles:

- Rational optimization of reprogramming culture conditions for the generation of induced pluripotent stem cells with ultra-high efficiency and fast kinetics
Jiekai Chen, Jing Liu, You Chen, Jiaqi Yang, Jing Chen, He Liu, Xiangjie Zhao, Kunlun Mo, Hong Song, Lin Guo, Shilong Chu, Deping Wang, Ke Ding, and Duanqing Pei
- The oncogene c-Jun impedes somatic cell reprogramming

Author information

Affiliations

Pei's Lab

Jing Liu , Jiekai chen & Duanqing Pei

Competing financial interests

The authors declare no competing financial interests.

Corresponding author

Correspondence to: Jing Liu (liu_jing@gibh.ac.cn)

Readers' Comments

Comments on this thread are vetted after posting.

Protocol Exchange ISSN 2043-0116

© 2015 Nature Publishing Group, a division of Macmillan Publishers Limited. All Rights Reserved.
partner of AGORA, HINARI, OARE, INASP, CrossRef and COUNTER

